

## Abbott Analytical



Consulting Scientists to the Disinfectant Industry

### Certificate of Analysis

Sample(s): One sample of Nilaqua

Received from: Waterless Limited

Hope Mill, Hope Mill business Centre, Brimscombe, Stroud,

GL5 2SE

Date received: 8 February 2010 Date tested: 11 February 2010

Certificate no: 10B.068ST4.CSS Certificate date: 19 February 2010 \*

**Sample ref:** 10B/068 **Page:** 1 of 2

Analysis required: Skin sensitivity test

#### Method:

Approximately 0.5 ml of Nilaqua was placed on a sterile lint gauze ( $2 \text{cm} \times 2 \text{cm}$  in size) and taped to the upper arm of each subject. A similar size piece of untreated gauze was also taped to the arm as a negative control. Each subject was asked to leave the gauze in place overnight and was examined the following day, comparing the treated patch with the untreated patch for any evidence of skin irritation.

The results are recorded in the table on page 2.

#### Conclusion:

Nilaqua does not show any significant adverse dermatological effect when in contact with human skin for  $24\ \mathrm{hours}$ .

D C Watson

<sup>\*</sup> Re-issued 17 February 2016 with amended product/brand names



# Abbott Analytical



### Consulting Scientists to the Disinfectant Industry

19 February 2010

Certificate No: 10B.068ST4.CSS

Page 2 of 2

#### Results:

Subject	Observations	
1	No redness	No irritation
2	No redness	No irritation
3	Slight redening of skin	No irritation
4	No redness	No irritation
5	No redness	No irritation
6	No redness	No irritation
7	No redness	No irritation
8	Slight reddening of skin	No irritation
9	No redness	No irritation
10	No redness	No irritation
11	No redness	No irritation
12	No redness	No irritation
13	No redness	No irritation
14	No redness	No irritation
15	Redness due to irritation	Slight irritation